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EXPLOITATION'S MANUAL

**PORATABLE APPARATUS FOR EHF- AND IR-THERAPY
WITH CHANGEABLE OSCILLATORS**

«AK TOM»

INTRODUCTION

Modern methods of treatment based on EHF- and BRR[®]-therapy in complex treatment helps to increase treatment effectiveness in many cases, as well as to decrease toxic load on the body. In some diseases, such as those concerned with the disturbance of blood microcirculation it helps to execute treatment in monotherapy mode.

EHF-therapy belongs to the class of non-invasive physiotherapy methods that influence the body in a «soft» and harmless manner, and increase its non-specific resistance. In some cases the «AK TOM» apparatus can be used in domestic conditions.

Modes of EHF-therapy are realized in the «AK TOM» apparatus with changeable oscillators. Hereinafter, under EHF-therapy the following types of radiation influence on the body are implied: background noise radiation and radiation at fixed frequencies.

Into the «AK TOM» oscillators realized the technology of background resonance radiation (BRR[®]), based on the application of materials that can reproduce own frequencies of the body in independent mode. Thanks to the influence on own frequencies of cell structures there is a weakening or a total suppression of pathological and intensification of physiological oscillations, with the restoration of dynamic balance in the patient's organism.

Extra high frequencies (EHF) occupy the spectrum of 30-300 GHz (wavelength spectrum is 10-1 mm). In the EHF-spectrum, the reaction of bio-objects is caused by frequency resonances, and high therapeutic effectiveness is recorded under weak radiation.

One of advantages of EHF-therapy include: painlessness (which is especially important for children and hypersensitive patients), the exclusion of accidental contamination of patients and medical staff.

The depth of penetration of EHF-radiation upon impact on the human skin surface is 300-500 μm , i.e. it is almost completely absorbed by the epidermis and top layers of the derma. The zone of its direct action covers nervous system receptors (mechanoreceptors, nociceptors and free nerve endings), the capillary channel of the blood system; skin leukocyte depot; cells of the diffuse neuroendocrine system, in particular mast cells. Therefore, the objects of therapeutic action of EHF-therapy are:

- nervous system;
- immune system and phagocytes;
- system for regulation of the blood's aggregative state (coagulative and anti-coagulative mechanisms);
- humoral control system.

As a result, the effectiveness of EHF-therapy is shown in a very broad spectrum of diseases.

Therapy with help of apparatus «AK TOM» is a new, state-of-the art method of EHF-therapy, which is based on technology for using materials with a controlled energy structure. They were created according to special technology on the basis of gallium arsenide, and capable of «remembering» external EHF-radiation, thus forming an exact similarity of signals of the pathogenic and/or medical factor. This enables an influence to be exerted on the organism not only in the classical mode of EHF-therapy, but also in the mode of background resonant radiation (BRR[®]). These are the modes executed by the «AK TOM» apparatus.

Clinical approbation has shown that EHF-therapy with help of BRR[®] is especially effective in cases of medication intolerance, presence of contraindications to methods of traditional physiotherapy and reflexotherapy, as well as when these forms of treatment prove insufficiently efficient.

The apparatus has passed clinical trials with numerous functional and laboratory studies that objectively registered a positive trend of treatment of many diseases.

The use of EHF-therapy in BRR[®] mode normalizes regulatory processes in the organism and it optimizes the rehabilitation of patients with acute and chronic diseases of various organs and organ systems, both at establishments for medicinal

and prophylactic treatment, and at home. Corresponding influence modes are specified additionally in «Physician's Compendium». The presented list does not limit indications for use of the apparatus, but each particular case should be preceded by consultation by an attending physician to determine the diagnosis and exclude individual contraindications.

ATTENTION!

- At the diseases that carry serious threats to health and life, the using of the device is held under the control of the doctor as a part of complex therapy.
- EHF-therapy has a gentle sedative effect, which, in event of individual heightened sensitivity to EHF radiation may lead to drowsiness and drop in attention span. In such an event work that requires concentration should be limited, including driving, after the procedure. Alternatively, the procedure should be applied before sleep.
- Is not allowed to use a mobile telephone or other means of wireless communication during the therapeutic procedure.
- Is not allowed to use the different types of oscillators at the same time.
- The total time of EHF-exposure in active mode per patient per day shall not exceed 30 minutes.

1. CONTRAINDICATIONS OF THE USE EHF-THERAPY

Contraindications of the use are:

- undetermined diagnosis;
- individual intolerance to this type of therapy;
- pregnancy;
- presence of implanted devices with autonomous power supply, such as a pacemaker and the like.

2. TECHNICAL DISCRIPTION OF THE «AK TOM» APPARATUS

«AK TOM» is the portable EHF-IR therapy apparatus with changeable oscillators. The oscillators are connected to the apparatus channels with the help of the flexible cables with sockets.

The apparatus «AK TOM» (Fig.1) has 11 modes of operations.

Each mode is designed taking into account effects on specific organs and body systems:

1 — «universal» EHF-mode;

2 — mode of «record» for the subsequent BRR®-influence;

- 3 — stressful conditions*;
- 4 — illnesses of the throat and larynx*;
- 5 — illnesses of the cardiovascular system*;
- 6 — illnesses of respiratory organs*;
- 7 — illnesses of endocrine system*, illnesses of digestive organs*;
- 8 — illnesses of musculoskeletal system and connective tissue*;
- 9 — illnesses of urinogenitals system*;
- 10 — illnesses of skin and subcutaneous fat*;
- 11 — operating mode on BAP*;


Fig. 1. Front view of the apparatus «AK TOM»

**- these modes are for use only by specialists or by the appointment of medical specialists.*

The presented list does not limit indications for use the apparatus, but each particular case should be preceded by consultation by an attending physician to de-

Contact for attachment of oscillator to the channel № 2

Contact for attachment of oscillator to the channel № 1

Display

Button for mode control

Button «Start»

Button to turne on the apparatus

The apparatus case


Fig. 2. Front view

Battery compartment
is closed

Battery compartment
is open

Fig. 3. Rear view

termine the final diagnosis and exclude individual contraindications. At the difficult diagnosis the initial course of treatment demands obligatory control and/or specialist participation.

The package contents of the apparatus comply with the Passport.

The apparatus «AK TOM» includes standard AA LR06 batteries.

Battery capacity is sufficient for approximately 3 month of everyday use.

IT IS CATEGORICALLY FORBIDDEN use salt batteries as the power supply elements (marked R06).

If the apparatus is not intended to be used for a long time, remove the batteries from the battery compartment.

If the battery charge has been exhausted, the icon appears on the display.

To change the batteries:

- open the battery compartment;
- take out the old batteries;
- insert new batteries into the apparatus;
- close the battery compartment.

Fig. 4. Changing the batteries

3. SETTING UP THE «AK TOM» APPARATUS

3.1. Insert two batteries into the apparatus (if they are not supplied with the apparatus).

3.2. Close the battery compartment.

3.3. Connect the oscillators to the apparatus using the cables.

3.4. Select the necessary mode by pressing the Mode button.

3.5. Image of the display in the selected mode shows information on the number of the mode and the length of the procedure in seconds.

3.6. Start the apparatus in the selected mode by pressing the START button.

3.7. When the apparatus is working in the selected mode, the display shows graphic and digital information on the duration of the procedure in second.

Evidence that the apparatus is working properly lies in the presence of marker (I and/or II, located in the lower line and which specifies that the channel is working) the changing time on the display and, simultaneously, the quiet «chirping» of the apparatus.

If the oscillator has failed or for any reasons it is not connected to the apparatus, the apparatus will «chirp», but the display will show nothing. After the procedure is complete, the «chirping» will stop, and the time indicator will arrive at zero.

Protect the apparatus and do not allow any liquid to get into the apparatus, especially inside the oscillators. The apparatus and oscillators can be disinfected only by wiping with cotton or a napking, slightly moistened with alcohol or a special disinfectant solution that causes no damage to plastic.

4. SEQUENCE OF ACTIONS WHEN WORKING WITH THE «AK TOM» APPARATUS

4.1. Connect the cable to the apparatus (to any channel). Connect the oscillator to the cable(Fig. 5).


Fig. 5

Both channels work independently from each other, therefore it is possible to connect to the apparatus as one or the both of oscillators at the same time (Fig. 6).


Fig. 6

4.2. Place the oscillator's working surface (smooth, white) on the skin's surface in the area of treatment (at open sores, wounds and burns may act through a bandage or the oscillator placed near the affected area). Can secure it with an adhesive bandage, leaving the cable free to be removed (Fig. 7).

Fig. 7



Fig. 8

4.3. Turn the apparatus on (Fig. 8).

4.4. Select the operating mode for the apparatus by pressing the mode buttons (Fig. 9).

Fig. 9

In this case on the indicator will be displayed:

- in the first line: the number of modes;
- in the second line: the time of the apparatus in this mode (in seconds).

4.5. Press the «Start» button to start the apparatus (Fig. 10).

In this case the unit will «chirr» and the display will show the following information:

- in the third line: timing scale;
- in the fourth line:
 - icon  if the oscillator connected to the first channel,
 - icon  if the oscillator connected to the second channel;
- in the second line will be shown the change of time (countdown showing how many seconds remain until the end of the operate in this mode).

Fig. 10

After a few seconds the indicator lights turn off, the apparatus will continue to work in this mode.

4.6. After the procedure is complete, the «chiring» will stop, the image on the display will return to its original state with the backlight off indicator. The apparatus goes into standby mode and after a few seconds automatically shut down.

To end a session early, simply press «Start» again.

While the apparatus is in standby mode, you can select a different mode by pressing the button.

The apparatus can work either in «Russian» or «Latin» mode. Switching from one mode to another is made in the «INFO» menu by choosing either «PYC» or «LAT» buttons (Fig. 11). After turning off the apparatus, the language chosen earlier will remain in place.


Fig. 11

Instead of adhesive bandage, you can use an elastic adhesive fasteners emitter, which does not irritate the skin and does not spoil the look of the oscillator.

Elastic mounting of the oscillator is supplied separately from the device.

Fig. 12

5. PRECAUTIONARY MEASURES AT WORK WITH THE APPARATUS «AK TOM»

5.1. You must be familiar with the current exploitation's manual before operating the apparatus.

5.2. The apparatus can be used in premises with an air temperature of between +10 °C and +35 °C and humidity of no more than 80 %.

5.3. While the apparatus is in operation, you must not expose unprotected eyes to the end of the oscillator that emits radiation (with the oscillator at a distance of less than 50 cm from the eyes). Doing so may result in undesirable physiological effects.

5.4. The apparatus is sensitive to the effects of electromagnetic discharge (ESD) on the connector cable to the detached oscillator. To prevent ESD, the following measures must be observed:

- connect and disconnect the oscillators only when the apparatus is turned off;
- before using the apparatus, the operator should discharge static electricity by touching a large metal object such as a heating pipe or water pipe;
- wear clothing that reduces the accumulation of electrostatic charge.

5.5. Before use, it is recommended that the oscillator be disinfected with a flannel cloth (or gauze, cotton ball), lightly moistened with ethyl alcohol. Before disinfection the cloth should be thoroughly squeezed out.

5.6. To prevent liquid from getting inside the body of the oscillator, **YOU MUST NOT** wipe the EHF oscillator with a very moist cloth

5.7. The «Mode» and «Start» buttons should only be pressed at intervals of at least 1—2 seconds, otherwise a failure or false start may occur.

5.8. Care must be taken when working with the detached EHF-oscillators to protect them from blows, bumps and falls to avoid damage.

5.9. The apparatus should only use the cable supplied, or a USB 2.0 AM/mini B 5P type cable from «Belsis», «Hama», «Sony», «Link Bits» or «Gembird». The length of the connector cable to the detached oscillator must be at least 0,5 m and no more than 1,8 m

5.10. The use of cables of other lengths or types is prohibited, as this can cause an increase in interference emission and a decrease in resistance to interference.

5.11. Connecting the apparatus and therapeutic oscillators to any other equipment via connecting cables is prohibited.

5.12. The oscillator connector cable must not be bent, broken, stretched while rubbing, or pulled when disconnecting the oscillator from the apparatus port.

5.13. When disconnecting or connecting the connecting cable to the apparatus or oscillator, the connector cable should be held only by the connector.

5.14. If the apparatus is not intended to be used for a long time, remove the batteries from the battery compartment.

To avoid damage to the apparatus and oscillators **YOU MUST NOT:**

- use salt batteries for the power supply (with the R06 marking);
- immerse the apparatus and the oscillators in water;
- subject the unit to blows and excessive mechanical stress;
- give the oscillators to children;
- swallow the oscillators.

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